



Oxivir® Tb

***Virucide, Bactericide, Tuberculocide,
Fungicide, Sanitizer**

ACTIVE INGREDIENT:

Hydrogen Peroxide 0.5%

OTHER INGREDIENTS: 99.5%

TOTAL: 100.0%

**KEEP OUT OF REACH
OF CHILDREN
CAUTION**

See additional precautionary
statements on side panel.



8 07174 49879 0



AHP® and Design are trademarks of Diversey, Inc.
USA Patent # 6346279, 6803057, 7632523,
8637085, 9233180

Net Contents

1 U.S. Qt. / 946 mL

SKU: 4277285



For Use in Health Care, Medical, Industrial & Institutional Facilities, Food Processing Establishments, and Hospitals.

This product cleans, disinfects and deodorizes hard, non-porous inanimate hospital or medical surfaces in one step with no rinsing required. Use daily on the following hard, non-porous inanimate surfaces: vinyl, painted surfaces, plastic surfaces, glazed ceramic, glazed porcelain, laminated surfaces and baked enamel surfaces associated with floors, walls, ceilings, tables, chairs*, non-food contact countertops, telephones, fixtures, glazed tile, toilet exteriors, urinal exteriors, sinks found in health care facilities, hospitals, patient rooms, operating rooms, shower rooms, and locker rooms. It can be used to pre-clean and disinfect hospital items: wheelchairs, patient bedrails and linings, patient bed linings, vinyl mattress covers, wash basins, and bed pans. A potable water rinse is required when disinfecting food contact surfaces. Do not use on glasses, dishes and utensils.

*Except Wood Surfaces

When used as directed, this product is effective against a wide variety of pathogenic microorganisms including bacteria, antibiotic resistant bacteria, viruses, and fungi.

BACTERICIDAL: (1 minute contact time) *Pseudomonas aeruginosa*, (ATCC 15442), *Staphylococcus aureus*, (ATCC 6538), *Salmonella enterica*, formerly known as *Salmonella choleraesuis* (ATCC 10708)

***VIRUCIDAL:** (1 minute contact time) *Norovirus, (Feline Calicivirus as the surrogate) (ATCC VR-782)

TUBERCULOCIDAL: in the presence of 5% serum load and 5 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces. *Mycobacterium bovis* (Tb) (OT 451C150)

FUNGICIDAL: (3 minute contact time) *Trichophyton interdigitale* (ATCC 9533) (formerly *Trichophyton mentagrophytes*), *Candida albicans* (ATCC 10231)

See reference sheet for a complete list of organisms and additional directions for use.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

For Use as a Cleaner or One-Step Disinfectant Product: Pre-clean heavily soiled areas. Spray evenly over hard, non-porous inanimate surface to be cleaned, 6 – 8 inches from surface, making sure to wet surfaces thoroughly. All surfaces must remain wet for 1 minute. Use a 5 minute contact time for Tb and a 3 minute contact time for fungi. Wipe surfaces dry, rinse or allow to air dry.

Note: Rinsing is not necessary unless floors are to be coated with finish or restorer. All food contact surfaces must be rinsed with potable water. Do not use on glassware, utensils, or dishes.

For Use as a Tuberculocide: Pre-clean heavily soiled areas. Spray evenly over hard, non-porous inanimate surface to be cleaned, 6 – 8 inches from surface, making sure to wet surfaces thoroughly. Allow the surface to remain visibly wet for 5 minutes. Wipe surfaces dry, rinse or allow to air dry. Change solution after each application.

To Kill Fungi: Pre-clean heavily soiled areas. Spray evenly over hard, non-porous inanimate surface to be cleaned, 6 – 8 inches from surface, making sure to wet surfaces thoroughly. Allow surface to remain visibly wet for 3 minutes. Wipe surfaces dry, rinse or allow to air dry.

To Spot Sanitize Soft Surfaces: Remove gross filth prior to sanitization. Hold bottle upright 6 – 8" from surface. Spray until fabric is wet. DO NOT SATURATE. Allow surface to remain visibly wet for 30 seconds. Gently blot with a clean cloth and let air dry. For difficult odors or heavy fabrics, repeat application. Repeat as needed.

Note: Prior to first use on a particular surface, it is recommended to test in a small, discreet area if there is concern of incompatibility.

For Use as a Cleaner and/or Deodorizer: Spray evenly over hard, non-porous inanimate surface to be cleaned, 6 – 8 inches from surface, making sure to wet surfaces thoroughly. Wipe surfaces dry, rinse or allow to air dry.

Manufactured for: © 2020 Diversey, Inc., PO Box 19747, Charlotte, NC 28219-0747 U.S.A. All Rights Reserved.

***KILLS HIV AND HCV and HBV ON PRE-CLEANED INANIMATE SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings (Hospitals, Nursing Homes) and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING AGAINST HIV-1 (AIDS Virus), HBV OR HCV OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable Protective Gloves (latex, nitrile, or vinyl gloves), gowns, face masks, and eye coverings must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/Body Fluids must be thoroughly cleaned from surfaces/objects before application of this product.

Contact Time: Allow surface to remain visibly wet for 1 minute to kill HIV-1, HBV and HCV. Use a 5 minute contact time for Tb and a 3 minute contact time for fungi.

Disposal of Infectious Material: Blood/Body Fluids must be autoclaved and disposed according to Federal, State and local regulations for infectious waste disposal.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal. Store in original container in areas inaccessible to small children. Keep securely closed.

PESTICIDE STORAGE: Store in a dry, well-ventilated area away from chemicals, direct light, heat or open flame. Do not mix with other cleaning or disinfecting products.

PESTICIDE DISPOSAL: Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility. For more information on the disposal of unused, unwanted product and the clean up of spills, contact the EPA Regional Office of the Manufacturer for guidance.

CONTAINER DISPOSAL: Nonrefillable container. Do not reuse or refill this container. Wrap empty container and put in trash.

ENVIRONMENTAL HAZARDS: This product is toxic to birds, fish and aquatic invertebrates. Caution should be used when applying indoors because pets may be at risk.

IN CASE OF EMERGENCY, CALL A POISON CONTROL CENTER OR DOCTOR FOR TREATMENT ADVICE.
1-800-851-7145

Have the product container or label with you when calling a Poison Control Center or doctor or going in for treatment.

EPA Reg. No. 70627-56

EPA Est. No. 74559-CAN-1 [VX]; 0312-WI-3 [NW]; 8251-OH-1 [KL]

Lot code letters indicate establishment number.

FOR INSTITUTIONAL USE

Questions? Comments:
1-800-558-2332

300203426 (20/163)

NOTE: Read and check carefully all information, including the specifications in this box. The job will not be sent to print until final okay is received. All information contained on the label and in this box are the responsibility of the final approver.




Patt. UV



PMS 2603



Black



PMS 185



PMS 072

Brand: Diversey

Product: Oxivir Tb

Pack Size: 1 Qt.

Date Started: 03/17/20

Label Size: 10.25 x 4.75

Printing: pressure sensitive label

Minimum Type: 6/6

UPC Code: UPC A

Product Code: 4277285

Tracking: 300203426

Replaces: 300203418

Coordinator: Corie Larson

Designer: M. Gary

Revision History:

03/24/20: Add EPA est number

03/30/20: Edits

05/07/20: Text edits

05/13/20: Edit Virucidal line

06/11/20: Remove errant "t"

SAP # 1280207