

OPTIM 33TB Wipes is an effective way to disinfect hard, non-porous environmental surfaces against:

***VRUCIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

* Poliovirus Type 1, Strain Chat

* Human Coronavirus

* SARS-Related Coronavirus 2 (SARS-CoV-2) (the virus that causes COVID-19)

* Influenza A virus, Strain A/8/68

* Herpes Simplex Virus, Type 1

* Herpes Simplex Virus, Type 2

* Rotavirus, Strain WA

* Rhinovirus, Type 14

* Kills HIV-1 [AIDS virus] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

* Kills HBV & HCV when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

* Kills Norovirus [Feline Calicivirus as the surrogate] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

This product has demonstrated effectiveness against SARS-CoV-2* (the virus that causes COVID-19).

HOSPITAL DISINFECTANT: BACTERICIDAL

in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Pseudomonas aeruginosa

Staphylococcus aureus

Escherichia coli O157:H7

Salmonella enterica

Acinetobacter baumannii

Staphylococcus aureus MRSA

Enterococcus faecalis VRE

TUBERCULOCIDAL in the presence of 5% serum load and 5 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Mycobacterium bovis

Non-Food Contact Surface Sanitizer for hard non-porous surfaces – in the presence of 5% serum load and 30 second contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Klebsiella pneumoniae

Salmonella enterica

Pseudomonas aeruginosa

Staphylococcus aureus

Enterococcus faecalis VRE

Escherichia coli O157:H7

Staphylococcus aureus MRSA

Notice to User: OPTIM 33TB Wipes may not be suitable for all items or surfaces. Prolonged exposure to copper, brass, anodized aluminum and some coatings should be avoided. If you are uncertain as to the material composition of your item, confirm with the manufacturer before proceeding.

For Use in Dental and Medical Offices

OPTIM 33TB Wipes is a one-step¹ disinfectant cleaner towelette designed for general cleaning and disinfecting of hard, non-porous environmental surfaces. OPTIM 33TB Wipes cleans quickly by removing dirt, grime, dead skin, blood, body soils and other organic matter commonly found in dental and medical offices.

OPTIM 33TB Wipes cleans, disinfects and deodorizes hard, non-porous environmental dental surfaces in one step¹. Its non-abrasive formula is designed for use on the following: painted surfaces, plastic surfaces, glazed ceramic, glazed porcelain, chrome, stainless steel, laminated surfaces and baked enamel surfaces associated with tables, chairs, countertops, telephones and sinks. OPTIM 33TB Wipes can be used to pre-clean and disinfect¹ non-critical medical equipment surfaces.

Directions for Use

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Not for personal cleansing. **This is not a baby wipe.** Do not flush in toilet.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

For Use as a Cleaner:

1. Pull towelette from canister and wipe hard, non-porous [environmental] surfaces.
2. Allow to air dry or rinse.

¹On hard, non-porous surfaces

²When used according to disinfection use directions.

³Use a 5 minute contact time for TB¹.

For Use as a One-Step Cleaner/Disinfectant:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute². Use additional towelettes if needed to ensure the correct wet contact time.
4. Allow to air dry or rinse.

For Use as a Non-Food Contact Sanitizer:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow surface to remain visibly wet for 30 seconds.
4. Allow to air dry or rinse.

NON UV AREA

AHP® and Design are trademarks of Diversey, Inc., used under license. OPTIM is a registered trademark of SciCan Ltd.

Distributed by:

SciCan Ltd., 1440 Don Mills Rd.
Toronto, ON M3B 3P9 Canada
Tel: 416 445 1600



Part No. **OPT33-W12** SCIL231 (25/038)

OPTIM 33TB WIPES

READY TO USE

One-Step¹ Disinfectant Cleaner

Tuberculocidal	5 MIN. AT 68°F (20°C)
*Virucidal	1 MIN. AT 68°F (20°C)
Bactericidal	1 MIN. AT 68°F (20°C)
Non-Food Contact Sanitizer	30 SEC. AT 68°F (20°C)

ACTIVE INGREDIENT

Hydrogen Peroxide	0.5%
Other Ingredients	99.5%
Total	100.0%

DO NOT FLUSH



KEEP OUT OF REACH OF CHILDREN

Net Contents:

1 pound (1lb.)
13 ounces (812 grams)
160 sheets 6" x 6.75"
(15.24 cm x 17.14 cm)

EPA Reg. No. 74559-3-83259
EPA Est. No. 74559-CAN-1

SciCan

STORAGE AND DISPOSAL:

Pesticide Storage: Store this product in a cool, dry area, away from direct sunlight and heat.

Container Handling: Dispose of wipe in trash after use. Do not flush. Tightly close lid flap between uses to retain moisture.

Non-refillable container. Do not reuse or refill this container. Offer empty container for recycling. If recycling is not available, discard container in trash.

Questions/Comments:

1-800-221-3046

Manufactured by:

Virox Technologies, Inc.
2770 Coventry Road Oakville, ON
L6H 6R1 Canada

*KILLS HIV-1 AND HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS

in health care settings and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/ OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, and eye coverings must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain wet for 1 minute² to kill HIV-1, HBV and HCV.

Disposal of Infectious Material: Cleaning materials used that may contain blood and other body fluids must be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

User should read the S.D.S. before using the product.
S.D.S. is available online at www.scican.com