



Titan™ Tabs

# Sporicidal Disinfectant Cleaner

## KEEP OUT OF REACH OF CHILDREN DANGER

**ACTIVE INGREDIENT:**

Sodium dichloro-s-triazinetriene ..... 48.21%\*

**OTHER INGREDIENTS:** ..... 51.79%

**TOTAL:** ..... 100.00%

\* Equivalent to 31.10% active chlorine by tablet weight. Refer to dilution chart for Available Chlorine concentrations.

<b>FIRST AID</b>	
<b>If in eyes</b>	<ul style="list-style-type: none"> <li>• Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li> <li>• Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.</li> <li>• Call a poison control center or doctor for treatment advice.</li> </ul>
<b>If swallowed</b>	<ul style="list-style-type: none"> <li>• Call a poison control center, or doctor immediately for treatment advice.</li> <li>• Have person sip a glass of water if able to swallow.</li> <li>• Do not induce vomiting unless told to do so by the poison control center or doctor.</li> <li>• Do not give anything by mouth to an unconscious person.</li> </ul>
<b>If inhaled</b>	<ul style="list-style-type: none"> <li>• Move person to fresh air.</li> <li>• If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.</li> <li>• Call a poison control center or doctor for further treatment advice.</li> </ul>
<b>If on skin or clothing</b>	<ul style="list-style-type: none"> <li>• Take off contaminated clothing.</li> <li>• Rinse skin immediately with plenty of water for 15-20 minutes.</li> <li>• Call a poison control center or doctor for treatment advice.</li> </ul>
<p><b>IN THE EVENT OF A MEDICAL EMERGENCY CALL YOUR POISON CONTROL CENTER AT 1-800-222-1222</b></p> <p>Have the product container or label with you when calling a poison control center or doctor, or going for treatment.</p>	
<p><b>NOTE TO PHYSICIAN</b></p> <p>Probable mucosal damage may contraindicate the use of gastric lavage.</p>	

See side panel for additional precautionary statements directions for use, storage and disposal and organisms.

For use in Cleaning and disinfection on hard, non-porous surfaces in schools, hotels, ice machines, spas, hot tubs, office buildings, hospitals, neo-natal units, hospital drain pipes, nursing homes, child care centers, daycares, veterinary clinics, kennels, boarding facilities, licensed care facilities, institutions, catering, vending machines, kitchens, Intensive Care Unit (ICU) operating rooms, dental facilities, gyms, health clubs, restrooms.

- Easy, Convenient, Accurate Measurement.
- Sustainable.

Convenient tablet form that is dissolved in water. Solution can be applied with dry wipes.

Can be applied with spray, mop, wipe, cloth, brush or mechanical sprayer.

Do Not Use with Hot Water!

**Net Contents:**

100 Tablets, 13.1 gram (0.46 oz.) Tablet Size

**SKU: 101102693**

Titan Tabs Sporicidal Disinfectant Cleaner provides effective cleaning strength that will not dull high gloss floors finishes with repeated use.

Penetrates biofilms, killing the bacteria\* living there.

Kills *C. diff* (*Clostridium difficile* spores) in 4 minutes!

Kills HIV-1, Hepatitis A Virus, Hepatitis B Virus and Hepatitis C Virus in 1 minute.

Can be used at a concentration of 1076ppm to disinfect floors, when used at the contact times specified in the Usage Table.

Cleans and Disinfects

### **Effervescent Disinfectant Tablets for Hospitals and Institutional Use:**

\*Titan Tabs Sporicidal Disinfectant Cleaner is effective against the following micro-organisms on pre-cleaned, hard, non-porous, inanimate surfaces: *Salmonella enterica*, *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Klebsiella pneumoniae*, *Staphylococcus epidermidis*, *Escherichia coli* O157:H7, *Staphylococcus aureus* – methicillin-resistant (MRSA) & glycopeptide-resistant (GRSA), carbapenem resistant *Klebsiella pneumoniae*, *Acinetobacter baumannii*, *Streptococcus pneumoniae*, vancomycin resistant *Enterococcus faecalis*, Poliovirus type 1, Herpes simplex virus type 1, Hepatitis A virus, Hepatitis B virus, Hepatitis C virus, Human Immunodeficiency Virus Type 1 (AIDS Virus), Influenza virus H1N1, respiratory syncytial virus, Canine Parvovirus, Newcastle Disease Virus, Pseudorabies, Canine Distemper Virus, Feline Calicivirus, Norovirus, Coxsackievirus, *Trichophyton interdigitale*, *Aspergillus fumigatus*, *Mycobacterium bovis* (TB) and *Clostridium difficile* spores. Refer to Usage Table for solution concentration and contact times.

## **PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

### **DANGER**

**Corrosive.** Causes irreversible eye damage. Harmful if swallowed, inhaled, or absorbed through skin. Do not get in eyes, on skin, or clothing. Avoid breathing dust. Wear chemical-resistant gloves and safety glasses or face shield when making up solution. Wash thoroughly with soap and water after handling, and before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove and wash contaminated clothing before reuse.

### **DIRECTIONS FOR USE**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Read the entire label and use strictly in accordance with precautionary statements and directions.

**General Solution Preparation:** Prepare a fresh solution twice weekly (every 3 days) with water of up to 400 ppm hardness when using closed containers (spray bottles). Prepare a fresh solution daily with water of up to 400 ppm hardness when using open containers (buckets) or if solution becomes diluted. Follow specific Directions for Use and Usage Table and Dilution Chart when preparing solution.

Do not use hot water in solution preparation.

All treated equipment that will contact food, feed, or drinking water must be rinsed with potable water before reuse.

**Notice to User:** This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

### Usage Table:

Pathogen	Minimum Dose required (ppm)	Minimum Contact time required (minutes)
<b>Sanitizer Claims</b>		
<i>Staphylococcus aureus</i> (ATCC 6538)	100 ppm	1 minute
<i>Salmonella enterica</i> (ATCC 6539)	100 ppm	1 minute
<b>Disinfection Claims – bacteria</b>		
<i>Staphylococcus aureus</i> (ATCC 6538)	538 ppm	10 minutes
	4306 ppm*	4 minutes
<i>Staphylococcus aureus</i> – methicillin resistant (MRSA) & glycopeptide-resistant (GRSA) (ATCC 33592)	1076 ppm	10 minutes
	4306 ppm*	4 minutes
<i>Staphylococcus epidermidis</i> (ATCC 51624)	1076 ppm	10 minutes
<i>Salmonella enterica</i> (ATCC 10708)	538 ppm	10 minutes
	4306 ppm*	4 minutes
<i>Pseudomonas aeruginosa</i> (ATCC 15442)	538 ppm	10 minutes
	4306 ppm*	4 minutes
<i>Streptococcus pneumoniae</i> (ATCC 6305)	1076 ppm	10 minutes
	4306 ppm*	4 minutes
<i>Escherichia coli</i> O157:H7 (ATCC 35150)	1076 ppm	10 minutes
<i>Acinetobacter baumannii</i> (ATCC BAA-1709)	4306 ppm*	4 minutes
Vancomycin resistant <i>Enterococcus faecalis</i> (ATCC 51575)	1076ppm*	10 minutes
	4306 ppm*	4 minutes
Carbapenem resistant <i>Klebsiella pneumoniae</i> (ATCC BAA-1705)	4306 ppm*	4 minutes
<i>Klebsiella pneumoniae</i> (ATCC 4352)	1076 ppm	10 minutes
<b>Biofilm Claims</b>		
<i>Pseudomonas aeruginosa</i> (in a biofilm) * (ATCC 15442)	4306 ppm	4 minutes
<i>Staphylococcus aureus</i> (in a biofilm) * (ATCC 6538)	4306 ppm	4 minutes

### Usage Table:

Pathogen	Minimum Dose required (ppm)	Minimum Contact time required (minutes)
<b>Virucidal Claims</b>		
Respiratory syncytial virus <sup>†</sup> (ATCC VR-26)	538 ppm*	10 minutes
Rhinovirus Type 14 <sup>†</sup> (ATCC VR-284)	1076 ppm*	10 minutes
Influenza Virus H1N1 <sup>†</sup> (ATCC VR-99) (ATCC VR-1469)	538 ppm* 4306 ppm*	10 minutes 1 minute
Human Immunodeficiency Virus Type 1 (HIV-1) <sup>†</sup> (Strain IIIB)	1076 ppm* 4306 ppm*	10 minutes 1 minute
Hepatitis A virus <sup>†</sup> (Strain HM175/18f)	1076 ppm 4306 ppm*	10 minutes 1 minute
Hepatitis B virus <sup>†</sup> (Duck Hepatitis B (DHBV)	1076 ppm* 4306 ppm*	10 minutes 1 minute
Hepatitis C virus <sup>†</sup> (Bovine Viral Diarrhea Virus Strain NADL – surrogate for Hepatitis C virus)	4306 ppm*	1 minute
Avian influenza A (H5N1) <sup>†</sup> (CDC #2006719965)	1076 ppm 4306 ppm*	10 minutes 1 minute
Norovirus <sup>†</sup> (ATCC VR-782)	2153 ppm*	1 minute
Poliovirus Type 1 <sup>†</sup> (ATCC VR-1000)	1076 ppm	10 minutes
Coxsackievirus B3 <sup>†</sup> (ATCC VR-30)	4306 ppm*	1 minute
Herpes simplex virus type 1 <sup>†</sup> (ATCC VR-733)	1076 ppm	10 minutes
<b>Fungicidal/Yeasticidal Claims</b>		
<i>Aspergillus fumigatus</i> (ATCC 36607)	4306 ppm*	1 minute
<i>Trichophyton interdigitale</i> (ATCC 9533)	1076 ppm	10 minutes
<b>Clostridium difficile Claims</b>		
<i>Clostridium difficile</i> spores (ATCC 43598)	2153 ppm 4306 ppm <sup>‡</sup>	10 minutes 4 minutes
<b>Mycobactericidal Claims</b>		
<i>Mycobacterium bovis</i> (TB) (ATCC 35743)	5382 ppm*	4 minutes

\*Note: testing has been conducted in the presence of 5% Serum soil load

‡Note: testing has been conducted in the presence of 0.25% Bovine Serum Albumin, 0.08% Bovine Mucin and 0.35% Yeast Extract soil load

## **HEALTHCARE and DISINFECTION/VIRUCIDAL† DIRECTIONS:**

Prepare solution strength as required, refer to Usage Table for correct doses and contact times; refer to Dilution Chart for solution preparation. Apply use solution to pre-cleaned, hard, non-porous, inanimate surfaces with mop, cloth, sponge, brush, wipe, or mechanical sprayer to wet all surfaces thoroughly. Allow surface to remain wet for contact time as indicated in the Usage Table, then remove product by wiping with brush, sponge, or cloth, or allow to air dry.

For sprayer applications using a spray device, spray at appropriate distance from surface depending on sprayer type (6 – 8 inches for spray bottles). Allow to remain wet for contact time as indicated in the Usage Table, then remove product by rubbing with brush, sponge, wipe or cloth or allow to air dry. Do not breathe spray mist.

Before using this product, food products and packaging materials must be removed from the room or carefully protected.

## **SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST Human Immunodeficiency Virus Type 1 (HIV-1) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:**

**PERSONAL PROTECTION:** Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks, and eye coverings.

**CLEANING PROCEDURE:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of Titan Tabs Sporidical Disinfectant Cleaner. This cleaning process may be accomplished with any cleaning solution including Titan Tabs Sporidical Disinfectant Cleaner.

**DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

**CONTACT TIME:** Leave surfaces wet for 10 minutes.

## **PERFORMANCE AGAINST BACTERIA GROWING IN A BIOFILM ON HARD NON-POROUS NON-FOOD CONTACT SURFACES**

Titan Tabs Sporidical Disinfectant Cleaner is also effective against bacteria ‡ *Staphylococcus aureus* and *Pseudomonas aeruginosa* growing in biofilms on hard, non-porous, non-food contact surfaces when used at a level of 4306 ppm available chlorine disinfectant solution with a 4 minute contact time.

## **DIRECTIONS FOR USE AGAINST BACTERIA GROWING IN A BIOFILM**

Pre-clean surfaces to remove soil and filth. Wipe dry. Prepare solution strength as required, refer to Usage Table for correct doses and contact times; refer to Dilution Chart for solution preparation. Thoroughly wet pre-cleaned surface with product. Allow surface to remain wet for 4 minutes. Rinse thoroughly.

## **DISINFECTION FOR SURFACES CONTAMINATED WITH CLOSTRIDIUM DIFFICILE**

Special Label Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores.

**Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks or eye covering.

**Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the disinfectant product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

**Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

## **Directions for Use:**

Prepare a 4,306 ppm solution; refer to Dilution Chart. Add 4 tablets to 1 gallon of water. Apply solution to pre-cleaned hard, non-porous, inanimate surface with mop, cloth, sponge, brush, wipe, or mechanical sprayer. Allow surface to remain wet for 4 minutes, then remove product by wiping with brush, sponge, or cloth, or allow to air dry.

## **HEALTHCARE DISINFECTION WITHOUT PRECLEAN**

Titan Tabs Sporidical Disinfectant Cleaner is a disinfectant when used at the doses and contact times indicated in the Usage Table.

## GENERAL and HEALTHCARE DISINFECTION WITHOUT PRECLEAN DIRECTIONS:

Prepare a 4,306 ppm solution; refer to Dilution Chart. Add 4 tablets to 1 gallon of water. Apply use solution to hard, non-porous, inanimate surface with mop, cloth, sponge, brush, wipe, or mechanical sprayer. Allow surface to remain wet for the appropriate contact time. Refer to Usage Table, then remove product by wiping with brush, sponge, or cloth, or allow to air dry.

## DISINFECTION FOR SURFACES CONTAMINATED WITH *MYCOBACTERIUM BOVIS* (TB) IN 4 MINUTES AT 20°C (68°F).

### Special Label Instructions for Cleaning Prior to Disinfection against *Mycobacterium bovis* (TB):

This product when used as directed below is effective against *Mycobacterium bovis* (TB) in 4 minutes at 20°C (68°F). This product can be used on hard non-porous surfaces in commercial institutional hospital and premises including kitchens, bathrooms, nurseries, sick rooms, laundry rooms, eating establishments, pet kennels, and veterinary premises. To disinfect hard non-porous surfaces, first clean surface by removing gross filth (loose dirt debris food materials etc). Prepare a 5,382 ppm available chlorine solution. Apply use solution to pre-cleaned, hard, non-porous, inanimate surfaces with mop, cloth, sponge, brush, wipe, or coarse mechanical sprayer to wet all surfaces thoroughly. Allow surface to remain wet for contact time as indicated in the Usage Table or for 4 minutes then remove product by wiping with brush, sponge, or cloth, or allow to air dry.

## SANITIZER PERFORMANCE

This product is an effective Sanitizer against *Staphylococcus aureus* and *Salmonella enterica* at 100 ppm with a 1 minute contact time.

## SANITIZING HARD, NON-POROUS SURFACES, DISHES, GLASSES, FOOD PROCESSING EQUIPMENT AND UTENSILS

Prepare a 100 ppm solution; refer to Dilution Chart for the number of tablets to use.

## STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

### PESTICIDE STORAGE

Store in a cool, dry, well-ventilated area at temperatures below 40°C/104°F. Avoid moisture getting into container

### PESTICIDE DISPOSAL

Pesticide may be acutely hazardous. Wastes resulting from the use of this product must be disposed of on-site, or at an approved waste disposal facility.

### CONTAINER HANDLING

Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available or puncture and dispose of in a sanitary landfill.

## PHYSICAL OR CHEMICAL HAZARDS

**STRONG OXIDIZING AGENT:** Use only clean dry utensils. Mix only into water. Contamination with moisture, dirt, organic matter, other chemicals or any other foreign matter may start a chemical reaction with generation of heat, liberation of hazardous gases and possible generation of fire and explosion. Avoid any contact with flaming or burning material such as a lighted cigarette. Do not use this product in any chlorinating device which has been used with any inorganic or unstabilized chlorinating compounds (e.g., calcium hypochlorite). Such use may cause fire or explosion.

## DILUTION CHART

Tablet size	13.1 g	
Solution ppm (mg/L) Available Chlorine	Tablets	Gallons of Water
0.5	1	2153
1	1	1076
1.5	1	718
3	1	359
4	1	269
5	1	215
10	1	100
100	1	10
538	1	2
1076	1	1
2153	2	1
4306	4	1
5382	5	1



**Questions? Comments**

**1-800-558-2332**

**EMERGENCY PHONE : 1-800-851-7145**

EPA Reg. No. 71847-7-70627

EPA Est. No. 71847-IRL-001



Manufactured for:  
Diversey Inc,  
PO Box 19747  
Charlotte, NC 28219-074

300205878 (19/149)



# Titan™ Tabs Sporidical Disinfectant Cleaner

**KEEP OUT OF REACH OF CHILDREN  
DANGER**

**DIRECTIONS FOR USE:**

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

See product label for additional precautionary statements, first aid, full directions for use and storage and disposal.

**Effervescent Disinfectant Tablets for Hospitals and Institutional Use:**

† Titan Tabs Sporidical Disinfectant Cleaner is effective against the following micro-organisms on pre-cleaned, hard, non-porous, inanimate surfaces: *Salmonella enterica*, *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Klebsiella pneumoniae*, *Staphylococcus epidermidis*, *Escherichia coli* O157:H7, *Staphylococcus aureus* – methicillin-resistant (MRSA) & glycopeptide-resistant (GRSA), carbapenem resistant *Klebsiella pneumoniae*, *Acinetobacter baumannii*, *Streptococcus pneumoniae*, vancomycin resistant *Enterococcus faecalis*, Poliovirus type 1, Herpes simplex virus type 1, Hepatitis A virus, Hepatitis B virus, Hepatitis C virus, Human Immunodeficiency Virus Type 1 (AIDS Virus), Influenza virus H1N1, respiratory syncytial virus, Canine Parvovirus, Newcastle Disease Virus, Pseudorabies, Canine Distemper Virus, Feline Calicivirus, Norovirus, Cocksackievirus, *Trichophyton interdigitale*, *Aspergillus fumigatus*, *Mycobacterium bovis* (TB) and *Clostridium difficile* spores. Refer to Usage Table for solution concentration and contact times.

**ACTIVE INGREDIENT:**  
Sodium dichloro-s-triazinetrione .....48.21%\*  
**OTHER INGREDIENTS:**.....51.79%  
**TOTAL:**.....100.00%  
\* Equivalent to 31.10% active chlorine by tablet weight. Refer to dilution chart for Available Chlorine concentrations

**Usage Table:**

Pathogen	Minimum Dose required (ppm)	Minimum Contact time required (minutes)
<b>Animal Pathogens<sup>1</sup></b>		
Canine Parvovirus* (ATCC VR-2017)	1076 ppm	10 minutes
Herpes simplex virus type 1* (ATCC VR-733)	1076 ppm	10 minutes
Newcastle Disease Virus* (ATCC VR-180)	1076 ppm	10 minutes
Pseudorabies* (ATCC VR-135)	1076 ppm	10 minutes
Feline Calicivirus* (ATCC VR-782)	1076 ppm	10 minutes
	2153 ppm*	1 minute
Canine Distemper virus* (ATCC VR-128)	1076 ppm*	10 minutes
Infectious Canine hepatitis* (ATCC VR293)	1076 ppm	10 minutes
Teschen/Talfan disease* (ATCC VR-669)	1076 ppm	10 minutes
Avian influenza virus (H5N1)* (ATCC VR-1608)	1076 ppm	10 minutes
	4306 ppm	1 minute
Porcine parvovirus* (ATCC VR-742)	1076 ppm	10 minutes
Runting & Stunting virus (tenosynovitis)* (ATCC VR-2449) ( ATCC VR-21)	1076 ppm	10 minutes
<i>Actinobacillus pleuropneumoniae</i> * (NCTC 12370) (ATCC 27088)	1076 ppm	10 minutes
<i>Bordetella bronchiseptica</i> (rhinitis)* (ATCC 19)	1076 ppm	10 minutes
<i>Brachyspira hyodysenteriae</i> (Treponema/Serpulina) (swine dysentery)* (ATCC 27164)	1076 ppm	10 minutes
Gumboro disease* (ATCC VR-478)	1076 ppm	10 minutes
<i>Streptococcus uberis</i> * (ATCC-9927)	1076 ppm	10 minutes
Transmissible gastroenteritis (TGE)* (ATCC VR-743)	1076 ppm	30 minutes
Swine Vesicular disease* (ATCC-VR-158)	1076 ppm	30 minutes
African swine fever* (ASFV)	1076 ppm	30 minutes
Hog cholera/Classical swine fever* (CSFV)	1076 ppm	30 minutes
Avipox (fowl pox)* (FPV)	1076 ppm	30 minutes
Respiratory syncytial virus* (ATCC VR-26)	538 ppm*	10 minutes
Bovine Viral Diarrhea Virus* (Strain NADL)	4306 ppm*	1 minute
Duck Hepatitis B Virus* (Duck Hepatitis B (DHBV)	4306 ppm*	1 minute
Porcine epidemic diarrhea virus* (Strain Colorado)	1076 ppm*	10 minutes

<sup>1</sup>Note: This use has not been approved by the California DPR  
<sup>2</sup>Note: these organisms not approved by the state of California  
<sup>3</sup>Note: testing has been conducted in the presence of soil load

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION IN ANIMAL HOUSING FACILITIES:**

1. Remove all animals and feed from premises, vehicles, and enclosures.
2. Remove all litter and manure from floors, walls and surfaces of barns, pens, stalls, chutes, and other facilities and fixtures occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding and watering appliances.
4. Thoroughly clean all surfaces with soap or detergent and rinse with water.
5. Saturate all surfaces with appropriate solution strength for the appropriate contact time, refer to Usage Table for correct dose and contact time, and to Dilution Chart for solution preparation
6. Immerse all halters, ropes, and other types of equipment used in handling and restraining animals, as well as forks, shovels, and scrapers used for removing litter and manure.
7. Ventilate buildings, cars, and other closed spaces. Do not house livestock or employ equipment until treatment has been absorbed, set, or dried.
8. Thoroughly scrub all treated feed racks, mangers, troughs, automatic feeders, fountains, and waterers with soap or detergent, and allow to air dry before reuse.

**EPA Reg No. 71847-7-70627**

**Establishment No. 71847-IRL-001**

Manufactured for: © 2019 Diversey, Inc., P.O. Box 19747, Charlotte, NC 28219-0747 U.S.A. All Rights Reserved.

REF90382 (19/149)