



ONE STEP • ONE MINUTE CLEANING & DISINFECTING

Virucidal • Bactericidal • Tuberculocidal • Fungicidal

ACTIVE INGREDIENTS

Ethyl Alcohol	72.50%
Didecyl Dimethyl Ammonium Chloride.....	.33%
OTHER INGREDIENTS*	27.17%
Total	100.00%

*Other ingredients do not include the weight of the dry wipes.

Net Contents:

**65 Pre-Moistened
Wipes (8" x 14")**
Net Wt: 1 lb. 10 oz. (726g)
Re-Order #: 60019542

**KEEP OUT OF REACH OF CHILDREN
CAUTION**

See side panels for additional precautionary statements

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals:

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

FIRST AID: If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after 5 minutes, then continue rinsing eye.

Call a Poison Control Center or doctor for treatment advice.

Have the product container or label with you when calling the Poison Control Center or doctor, or going for treatment. You may also contact Intotrac 24-Hour Emergency Service at 1-800-858-7378 seven days a week 6:30 am to 4:30 pm Pacific Time (NPIC Web Site www.npic.orst.edu)

PHYSICAL OR CHEMICAL HAZARDS: Flammable. Keep away from heat or open flame.

STORAGE AND DISPOSAL: Do not contaminate water, food or feed by storage or disposal.

PESTICIDE STORAGE: Always store this product and all pesticides in the original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold.

PESTICIDE DISPOSAL: Waste resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

CONTAINER DISPOSAL: Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in trash.

Ultracare® Minit Wipes are intended for use in:

- Hospitals and healthcare settings, including: Ambulances, Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's Office, Emergency Rooms, Food Service and Food Prep Areas, Home Health Care, Hospices, Laboratory, Laundry Rooms, Long Term Care Facilities, Newborn Nursery, Nursing Homes, Offices, Operating rooms, Pediatrics, Physical Therapy, Physician's Office, Ophthalmic Offices, Radiology, Recovery Room, Rehabilitation, Respiratory Therapy, Ultrasound, Transport Vehicles, Waiting Areas, X-Ray.
- Critical Care Areas: (CCU) Cardiac Care Units, (E.R.) Emergency Rooms, (ICU) In-

tensive Care Units, (NICU) Neonatal Intensive Care Units, (PICU) Pediatric Intensive Care Units, Surgery, Surgical Prep, and Surgical Recovery Units.

- Animal Holding Areas, Animal Hospitals, Animal Life Science Laboratories, Animal Transportation Vehicles, Breeding Establishments, Grooming Establishments, Kennels, Pet Shops, Tack Shops, Veterinary Clinics, Veterinary Examination Rooms, Veterinary Operating Rooms, Veterinary Waiting Rooms, Veterinary X-Ray Rooms, Zoos.
- Childcare Facilities, Daycare Centers, Nurseries.
- Households, Homes, Residences, Kitchens and/or Bathrooms.
- Institutional Kitchens and/or Bathrooms.

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not for use on critical and semi-critical medical device surfaces.

CLEANING: Use one or more wipes, as necessary, to wet surface sufficiently and to thoroughly clean the surface.

DISINFECTING: To disinfect hard, non-porous surfaces, use one or more wipes, as necessary, to thoroughly wet the surface to be treated. Treated surface must remain visibly wet for one (1) minute to achieve complete disinfection (or kill) of all pathogens listed on this label. For visibly soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry. For surfaces that may come into contact with food, a potable water rinse is required.

***VIRUCIDAL:** Ultracare® Minit Wipes are an effective virucide on hard, non-porous surfaces, in the presence of an organic load, against Poliovirus serotype 1 (PV-1) [ATCC VR-1562], Hepatitis B Virus, Hepatitis C Virus, Rhinovirus [ATCC VR-340], Norovirus, Human Coronavirus [ATCC VR-1558], Rotavirus Strain WA, Respiratory Syncytial Virus [ATCC VR-29], Human Immunodeficiency Virus (HIV-1), Herpes Simplex Virus Type 1 [ATCC VR-539], and Influenza Virus type A2 [ATCC VR-544], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

BACTERICIDAL: Ultracare® Minit Wipes are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against Escherichia coli (E. coli) O157:H7 [ATCC 35150], Pseudomonas aeruginosa [ATCC 15442], Salmonella enterica (Salmonella) [ATCC 10708], Staphylococcus aureus [ATCC 6538], Streptococcus pyogenes (Strep)

[ATCC 19615], Vancomycin resistant *Enterococcus faecium* [VRE] [ATCC 51559], Methicillin resistant *Staphylococcus aureus* [MRSA] [ATCC 33591], *Klebsiella pneumoniae* [ATCC 4352], *Listeria monocytogenes* [ATCC 984], *Acinetobacter baumannii* [ATCC 19606], *Campylobacter jejuni* [ATCC 29428], and *Burkholderia cepacia* [ATCC 25416], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

TUBERCULOCIDAL: Ultracare® Minit Wipes are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Mycobacterium bovis* [ATCC 35743], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

FUNGICIDAL: Ultracare® Minit Wipes are an effective fungicide against *Trichophyton mentagrophytes* [ATCC 9533], causative agent of Athlete's Foot Fungus, when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for 1 minute. Allow surfaces to air dry.

KILLS HIV-1, HBV and HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV and HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye coverings.

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of Ultracare® Minit Wipes.

CONTACT TIME: Allow surface to remain wet for 1 full minute.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

EFFICACY TESTS HAVE DEMONSTRATED THAT ULTRACARE® MINIT WIPES ARE AN EFFECTIVE, VIRUCIDE, BACTERICIDE, TUBERCULOCIDE AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL.

DISPENSING INSTRUCTIONS:

To Dispense Wipes: Open lid, remove entire cover, remove and discard seal. From center of wipes roll pull up a corner of the wipe, twist it and thread it through the hole in the container cover. Pull through about one inch. Replace cover tightly. Pull out first wipe and tear off at a 45° angle. Remaining wipes feed automatically, ready for next use. When not in use, keep the tip of the exposed wipe inside the well in the cover, and the lid snapped closed.

No warranty, either expressed or implied, is given to the usage contrary to the instructions given herein.

Veterinary Clinics, Animal Life Science Laboratories, Pet Shops : Remove all animals and feed from premises, vehicles, and enclosures. Remove all litter and manure from floors, walls, and surfaces of facilities and fixtures occupied or traversed by animals. Empty all troughs, racks, and other feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Use one or more wipes to wet surfaces for 1 minute. Immerse all handling restraining equipment as well as forks, shovels, and scrapers used for removing litter and manure. Ventilate buildings, vehicles, and other closed spaces. Do not house animals or employ equipment until treatment has been absorbed, set, or dried. Thoroughly scrub treated feed troughs, racks and other feeding and watering appliances with soap or detergent, then rinse with potable water before reuse.

Ultracare® Minit Wipes will not harm most surfaces: including stainless steel, aluminum, plastics, fiberglass, vinyl, and most painted surfaces. For plastic and painted surfaces, spot test on an inconspicuous area before use. Do not use on unpainted wood. Not for use on people or pets.

EPA Reg. No. 88494-2-63562
EPA Est. No. 85743-TUR-1
MADE IN Turkey

Manufactured for:
Dentronix
235 Ascot Parkway | Cuyahoga Falls, OH 44223 | U.S.A.
800.523.5944 | Fax 1.330.916.7077
info@dentronix.com | www.dentronix.com

Ultracare®
By: dentronix®



This container is coded for recycling.



+D661600195420

P/N 31233A